

# POLYNUCLEOTIDE PHOSPHORYLASE (PNPase)

[EC 2. 7. 7. 8]

from Bacillus stearothermophilus

 $RNA_{n+1} + Pi \leftrightarrow RNA_n + Nucleoside diphosphate$ 

# FOR DEPOLYMERIZATION REACTION

SPECIFICATION		
State	: Lyophilized	
Specific activity	: more than 2,000 U/mg protein	
PROPERTIES		
Molecular weight	: 300,000 - 340,000	
Subunit molecular weight	: ca. 85,000	
Optimum pH	: 9.0 - 9.5	(Fig. 1)
pH stability	: 9.0 - 11.0	(Fig. 2)
Isoelectric point	: 4.0	( <b>U</b> )
Thermal stability	: No detectable decrease in activity up to 55 °C.	(Fig. 3, 4)
Michaelis constants	: (38 mM Tris-HCI buffer, pH 9.5, at 60 °C)	
	Poly A	0.27 mM**
	KH <sub>2</sub> PO₄	3.0 mM
	**concentration of poly A was calculated as AMF	concentration
Effectors	: cations and anions	(Fig. 5, 6)

## STORAGE

Stable at -20 °C for at least one year

# APPLICATION

The enzyme is useful for the preparation of polyribonucleotide.



## ASSAY

# Principle

The change in absorbance is measured at 340 nm according to the following reactions.

Poly  $A_n + Pi$  <u>PNPase</u> Poly  $A_{n-1} + ADP$  (I) ADP + PEP <u>PK</u> ATP + PyruvatePyruvate + NADH + H<sup>+</sup> <u>LDH</u> Lactate + NAD<sup>+</sup> (II)

# **Unit Definition**

One unit of activity is defined as the amount of PNPase that forms 1  $\mu$ mol of ADP per hour at 60 °C by depolymerizing of Poly A.

# Solutions

(Reaction I)

- I Buffer solution ; 100 mM Tris-HCl, pH 9.5 ((1.212 g Tris + 0.074 g EDTA + 0.014 mL 2-mercaptoethanol + 0.610 g MgCl<sub>2</sub>· $6H_2O$  + 0.746 g KCl)/80 mL distilled water, adjusted to pH 9.5 with 1 N-HCl and filled up to 100 mL with distilled water)
- II KH<sub>2</sub>PO<sub>4</sub> solution ; 65 mM (0.088 g KH<sub>2</sub>PO<sub>4</sub>/10 mL distilled water)
- III polyadenylate (Poly A) solution ; (25 mg Poly A potassium salt/1 mL distilled water; ca. 35 mM based on AMP concentration)

# (Reaction II)

- IV Buffer solution ; 100 mM Triethanolamine buffer, pH 7.6 ((9.300 g triethanolamine-HCI + 0.407 g MgCl<sub>2</sub>·6H<sub>2</sub>O + 0.373 g KCI)/400 mL distilled water, adjusted to pH 7.6 with 1 N-NaOH and filled up to 500 mL with distilled water)
- V NADH solution ; 13.1 mM (0.100 g NADH disodium salt·3H<sub>2</sub>O/10 mL distilled water)
- VI Phosphoenolpyruvate (PEP) solution ; 56mM (0.150 g PEP MCA salt/10 mL distilled water)
- Ⅶ Pyruvate kinase (PK); (from rabbit muscle, Roche Diagnostics K.K., No. 128 155) crystalline suspension in 3.2 M (NH₄)<sub>2</sub>SO₄ solution (10 mg/mL) approx. 200 U/mg at 25 °C
- Ⅷ Lactate dehydrogenase (LDH); (from pig heart, Oriental Yeast Co. Ltd., Product Code: LDH-02) ammonium sulfate suspension, approx. 5,000 U/mL at 25 °C

## Preparation of Enzyme Solution

Dissolve the lyophilized enzyme with distilled water and dilute to 1 to 5 U/mL with 50 mM Tris-HCl buffer, pH 8.5.

## Procedure

(Reaction I)

1. Prepare the following reaction mixture and pipette 0.55 mL of reaction mixture into a test tube.

Solution I	2.50mL	Solution III	1.00mL
Solution II	1.00mL	H <sub>2</sub> O	1.00mL

- 2. Add 0.10 mL of enzyme solution and mix.
- 3. Incubate at 60 °C for exactly 10 minutes.
- 4. After incubation, add 0.01 mL conc. HCI and mix.
- Centrifuge at 10,000 rpm for 30 seconds. At the same time, repeat the Procedure 1 to 5 using distilled water in place of enzyme solution in Procedure 2 (as blank).

(Reaction II)

6. Prepare the following reaction mixture and pipette 2.50 mL of the reaction mixture into a cuvette. Solution IV 24.18mL Solution VII 0.12mL



Solution V	0.40mL	Solution	0.05mL
Solution VI	0.25mL		

- 7. Incubate at 30 °C for about 3 minutes.
- 8. Add 0.10 mL of supernatant of Procedure 5 and mix.
- 9. Read absorbance at 340 nm (Abs•test). Repeat the Procedure using blank (Abs•blank).

## Calculation

Volume activity (U/mL) = ((Abs•blank) - (Abs•test)) X  $\frac{2.60 \times 0.65}{6.22 \times 0.10 \times 0.10} X \frac{60}{10} X d.f.$ 

Volume activity (U/mL)

Specific activity (U/mg protein) =

Protein concentration (mg/mL)\*

d.f.; dilution factor

6.22 ; millimolar extinction coefficient of NADH ( $cm^2/\mu mol$ ) \*Protein concentration ; determined by the absorbance at 280nm (Abs280), where 1 Abs280 = 1 mg/mL

## REFERENCES

- 1. Smith, J.C., and Eaton, M.A.W.; Nucleic Acids Research, 1, 1763 (1974)
- 2. Wood, J.N., and Hutchinson, D.W.; ibid., 3, 219 (1976)









JIPZ

NIPRO ENZYMES

Measurement : 0.015 mL of each cation solution, 0.010 mL of enzyme solution and 0.055 mL of reaction mixture were mixed, and reacted at 60 ° C. After 10 minutes, the quantity of ADP was determined. O NaCl,  $\triangle$  KCl,  $\Box$  MgCl<sub>2</sub>,  $\blacksquare$  CaCl<sub>2</sub>,  $\blacktriangle$ ZnCl<sub>2</sub>



Fig. 6 Effect of various anions on the activity of Polynucleotide phosphorylase in the following Assay Method

Measurement : 0.015 mL of each anion solution, 0.010 mL of enzyme solution and 0.055 mL of reaction mixture were mixed, and reacted at 60 ° C. After 10 minutes, the quantity of ADP was determined. O NaCl,  $\triangle CH_3COONa$ ,  $\Box Na_2SO_4$ ,  $\blacksquare NaHCO_3$ ,  $\blacktriangle NaH_2PO_4$