

# PHOSPHOTRANSACETYLASE (PTA)

# [EC 2. 3. 1. 8]

from Bacillus stearothermophilus

Acetyl-CoA + Pi ↔ Acetylphosphate + CoA

SPECIFICATION		
State	: Lyophilized	
Specific activity	:more than 5,000 U/mg protein	
Contaminants	:(as PTA activity = 100 %)	
	Acetate kinase	< 0.01 %
	Adenylate kinase	< 0.01 %
	Lactate dehydrogenase	< 0.01 %
PROPERTIES		
Molecular weight	: ca. 70,000	
Subunit molecular weight	: ca. 35,000	
Optimum pH	: 7.5	(Fig. 1)
pH stability	: 7.0 - 11.0	(Fig. 2)
Isoelectric point	: 4.5	
Thermal stability	: No detectable decrease in activity up to 50 °C.	(Fig. 3, 4)
Michaelis constants	:(87 mM Tris-HCl buffer, pH 7.5, at 30 °C)	
	Coenzyme A	0.4 mM
	Acetyl Phosphate	1.1 mM

# STORAGE

Stable at -20 °C for at least one year

#### APPLICATION

The enzyme is useful for determination of CoA or acetate.



# ASSAY

#### Principle

The change in absorbance is measured at 233 nm according to the following reaction.

Acetylphosphate + CoA \_\_\_\_\_ Acetyl-CoA + Pi

#### **Unit Definition**

One unit of activity is defined as the amount of PTA that forms 1  $\mu mol$  of acetyl-CoA per minute at 30  $^\circ C.$ 

#### Solutions

- I Buffer solution ; 100 mM Tris-HCl, pH 7.5
- I CoA solution ; 6.4 mM (50 mg CoA trilithium salt/10 mL distilled water)
- III Acetylphosphate solution ; 217 mM (0.400 g acetylphosphate potassium lithium salt/10 mL distilled water)
- IV Ammonium sulfate (AmS) solution ; 1 M (13.2 g AmS/100 mL distilled water)

# **Preparation of Enzyme Solution**

Dissolve the lyophilized enzyme with distilled water and dilute to 5 to 20 U/mL with 50 mM Tris-HCI buffer, pH 8.0.

#### Procedure

- 1. Prepare the following reaction mixture and pipette 3.00 mL of reaction mixture into a cuvette.
  - Solution I 26.0 mL Solution II 1.0 mL
    - Solution II 2.0 mL Solution IV 1.0 mL
- 2. Incubate at 30 °C for about 3 minutes.
- 3. Add 0.01 mL of enzyme solution into the cuvette and mix.
- 4. Read absorbance change at 233 nm per minute ( $\Delta Abs_{233}$ ) in the linear portion of curve.

# Calculation

Volume activity (U/mL) = 
$$\frac{(\Delta Abs_{233}) X (3.00 + 0.01)}{4.44 X 0.01} X d.f.$$

Specific activity (U/mg protein) = Volume activity (U/mL)

Protein concentration (mg/mL)\*

d.f. ; dilution factor

4.44 ; differential millimolar extinction coefficient between acety-CoA and CoA (cm<sup>2</sup>/µmol) \*Protein concentration ; determined by Bradford's method



