DENTAL NEEDLE

Instructions for Use

Prior to use, please read the Instructions for Use carefully.


[Indications]

The device is used by a dentist to anesthetize the site of a patient that requires painful dental treatment.

[Contraindications]

- Do not re-use

[Shape and Structure]

[Clinical benefits]

By connecting to a dental cartridge syringe, an anesthetic can be administered to the gums, which act as a travel pathway for peripheral nerves, thus reducing the pain experienced by patients who undergo dental treatment.

[Cautions]

1. The needle tube may break if the needle is bent repeatedly during use, so use with caution.
2. Stop using the device if damage such as a hooked needle tip occurs.
3. The needle might cut the rubber of the cartridge when it is connected to the cartridge syringe. Perform a test injection or other procedure to test for clogging of the needle.
4. Handle the device carefully so as to not stick yourself unintentionally.
5. Do not apply excessive force to the needle hub. Also, do not apply excessive rotational force through the container when connecting the needle hub. The needle tube or needle hub could break.
6. When removing the container or container cap, ensure that they do not come into contact with the needle tip. The needle tip could deform and make puncturing difficult or block the flow. Pull straight out.
7. Do not allow drugs to get on the connecting parts. They may loosen.
8. Do not recap the needle. Recapping may result in needle sticking. Also, if the container is incorrectly recapped at an angle, the container or container cap may be pierced by the needle.
9. This is a synthetic resin device that deteriorates in low temperatures. Handle carefully in low temperature environments. It may break.
10. Re-use of the device may lead to infection by contamination.
11. Dispose of the device in an approved biohazard container as per facility protocol.
12. Any incidents shall be reported to the manufacturer and the competent authority of your State.
13. Do not use if package is damaged.
14. If the device is exposed to abnormal conditions (for example, high temperature and humidity) or if it is unintentionally opened before use, do not use the device.

[Instructions for Use]

1. Twist the container to break the sealing paper, and then remove the container cap.
2. Screw the needle hub onto the dental cartridge syringe to secure it in place.
3. Remove the container, and proceed with using the device.

[Connection with other devices]

The device connects to a dental cartridge syringe with the screw part of the needle hub.

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