SUREFLUX™ Series FB Series

HEMODIALYZER

Instructions for Use

Indications: This device is indicated for patients having acute or chronic renal failure when dialysis is prescribed by the physician.

Contraindications: Do not re-use. Patients indicating allergic reactions to cellulose triacetate membranes should not be dialyzed with this product.

Clinical Benefits: Dialysis treatment can be performed by connecting the device to a dialysis machine, a blood line, and an AVF needle, and waste and unnecessary water is removed from the blood to support the life of renal failure patients.

Warning

- Use this product according to instructions of a physician who is well familiar with the patient's condition.
- Follow these instructions and those of the dialysis machine supplier.
- Do not use for any other purposes than dialysis.
- If any abnormalities such as foam generation or mixture, blood leakage, blood coagulation and hemolysis occurred during the use of this product, take appropriate measures according to a physician's instructions.
- In case drugs including an anticoagulant are administered before or during use of this product, follow a physician's instructions about the administration and dose and the administration time of the drugs.
- Do not re-use this product since this is a single-use product.
- The foreseeable risks in association with re-use of the product are:
 - Infection by contamination,
 - Deterioration of solute removal performance and ultrafiltration performance,
 - Exposure of patients and/or technicians to residual medicinal agents such as disinfectant used for product re-use, and/or adverse effects of residual medicinal agents on them, and
 - Damage of hollow fiber and/or leakage.
- Do not expose this product to chemical solvents, such as bleach and alcohols. Residual disinfectant in this product may cause adverse patient reactions.
- If the patient exhibits any abnormal symptoms such as discomfort, pruritus, urticaria, peripheral and facial edema, respiratory arrest, facial flush, erythema, asthmatic reaction, hypertension, hypotension and/or cardiac arrhythmia during the use of this product, take appropriate measures according to a physician's instructions.
- Commonly seen side effects (hypotension, hypertension, headache and nausea which are sometimes with hypovolemia or hypervolemia) can be minimized by careful management of the patient fluid and electrolyte balance, as well as the dialysis condition (blood flow rate and ultrafiltration rate).
- During dialysis, constantly monitor the patients who:
 - have a history of hypotension with therapy,
 - have inflammatory reaction, allergic reaction, hypersensitivity, or increase in the immunity by infections,
 - take hypotensive drugs such as inhibitor of angiotensin converting enzyme and calcium antagonist,
 - use this product for the first time.

Caution: Caution should be employed against excessive water removal. Use of accurate UF control system is required. Do not use on non-de-aerated dialysis fluid delivery systems. Confirm that the dialysate does not contain pyrogens in order to prevent transfer of pyrogens from dialysate to blood.

Any incidents shall be reported to the manufacturer and the competent authority of your State.

Caution before use

- Do not use if package is damaged.
- Do not use the device if the protective end caps are not in place.
- Unpack immediately before use.
- Avoid air mix-in and contamination during rinsing / priming operations.
- Start dialysis immediately after rinsing / priming operations.
- Rinsing / priming should be carried out under the following conditions according to this "Instructions for Use":
 - Blood side : Rinse with physiological saline (not less than 1 000 mL) at a flow rate of 100 mL/min.
 - Dialysate side : Verify conductivity and temperature, and rinse with dialysate at a flow rate of 500 mL/min for about 5 minutes.
- Check the integrity of the blood line and dialyzer.
- Administration of Heparin Systemic or regionalized heparinization may need to be administered based on instructions from the attending physician.

Caution in use

- Continuously monitor the pressure in the blood line and check for blood leakage during dialysis.
- Carefully avoid contamination during blood sampling and blood reinfusion.
- Set TMP alarm (max. 500 mmHg).
- Avoid air during blood reinfusion to minimize the risk of air embolisms.
- Do not apply excessive pressure to the blood line, the dialyzer and their connections.

Caution after use

- Single use only. Dispose of the dialyzer immediately after use.
- Dispose of the used blood lines and dialyzer by any means suitable for avoiding contamination.
- Dispose of the device in an approved biohazard container as per facility protocol.

Caution for storage

- If the device is exposed to abnormal conditions (for example, high temperature and humidity) or if it is unintentionally opened before use, do not use the device.
- Handle carefully in low temperature environment. (Especially below freezing, it may result in breakage of the hollow fiber.)

Instructions for Use

I. Rinsing / priming

Follow instructions on the machine's operator manual.

«In case of priming with a physiological saline bag»

- Take the dialyzer out of the package and set it to the holder so that the venous side is directed upward. (Fig. 1)
- Connect the arterial and venous dialyzer connector to dialyzer. Connect the arterial patient line to physiological saline bag (Fig. 2). Start arterial blood pump at a flow rate of 100 mL/min (not less than 1 000 mL) (Fig. 3).
- Stop the blood pump; rotate the dialyzer for 180 degree. Put the dialysate connectors to the dialyzer (dialysate inlet at the venous blood side, dialysate outlet at the arterial blood side) (Fig. 4). Make sure of degassing the dialysate part of the dialyzer. Start running through dialysate at a flow rate of about 500 mL/min.
- Restart the blood pump. Make sure that the blood compartment is free of air-bubbles and filled with physiological saline. The preparations for dialysis is complete (Fig. 5).

«In case of online priming»

Follow instructions on the machine's operator manual.

[Leakage test] In case of physiological saline priming

It is recommended to perform the following operations before connecting the dialysate lines to the dialyzer.

- Fully prime the arterial and venous lines and the dialyzer with physiological saline by operating the blood pump ; then stop the pump operation.
- Clamp the arterial line near the dialyzer and the distal end of the venous line with clamps.

- Place the clamped distal end about 1 m below the dialyzer and remove the clamps. (This results in application of a negative pressure of approximately 70 mmHg to the blood compartment of the dialyzer.)

- Examine whether or not continuous bubble formation is observed in the venous header to check for leakage from the dialyzer, if observed, replace the dialyzer with a new one.

II. Start of dialysis

Follow instructions on the machine's operator manual.

- Prepare the blood access site and connect to the arterial line. Remove the clamps from the arterial and venous lines. While running dialysate at a flow rate of approximately 500 mL/min, operate the blood pump at a flow rate of approximately 50 mL/min.
- Confirm that no air bubbles remain in venous header or venous blood line.
- Fully prime the arterial and venous lines including the dialyzer with blood by operating the blood pump, and then stop the pump operation. Clamp the distal end of the venous line with clamps.
- Prepare the blood return site and connect to the venous line. After confirming that there are no bubbles in the line, remove the clamps from the line. After checking that there are no clamps on the lines and no line folding, operate the blood pump at a low flow rate. Check the integrity of the connections.
- After confirming that there are no bubbles in the arterial and venous headers, turn the dialyzer 180° to allow removal of bubbles from dialysate. If bubbles are detected in the venous header before the turning, run blood at a prescribed flow rate for 5 to 10 minutes with the venous side kept upward.

III. Operations during dialysis

- If stopping the blood pump is required during dialysis due to insufficient blood flow, lower the dialysate pressure to approximately 0 mmHg. (This is to avoid blood coagulation due to dehydration.)
- Set the UF rate carefully to avoid excessive water removal according to the patients' needs. Reduce the blood flow rate if disequilibrium syndrome is suspected.
- If blood leakage is suspected, judge by testing dialysate sampled from the dialysate outlet port of the dialyzer using occult blood reaction test paper. If a leakage is detected, reduce the UF rate to minimum rate according to the institutional protocol, stop the dialysate supply and blood reinfusion, then replace the dialyzer with a new one.

IV. Dialysis termination and blood reinfusion

- Stop the blood pump, clamp the arterial line and remove the line from the arterial blood access site ; then connect the line to the physiological saline vial for blood reinfusion.
- Unclamp the arterial line and run physiological saline to rinse out the blood from the arterial and venous lines and the dialyzer.
- After the blood reinfusion, discard the arterial and venous lines and the dialyzer. Do not reuse them.

Performance and Specification

The performance of the hemodialyzer varies with types. Refer to respective catalogues and performance data sheets.

Connection with other devices

The device contacts and is connected with the blood line at the header and the dialysis machine through the coupler port of the housing.

Guarantee

- Non-pyrogenic, sterile fluid path that has been sterilized using irradiation.
- This dialyzer is manufactured under strict quality control and the quality is assured. If the dialyzer is defective (broken package, damaged dialyzer), it shall be replaced with a new one at our cost upon return of the broken package or damaged dialyzer. We will not be responsible, however, for the injury of a patient or any person or the damage to any object that is attributed to transport, storage and operation in your institution.
- If a patient or any person is injured or any object is damaged by use of our dialyzer, we will not be responsible for the injury or damage unless we are clearly identified as being at fault.
- If a patient or any person is injured or any object is damaged by reuse of our dialyzer, we will not be responsible for the injury or damage of any nature.
- We will not be responsible for any injury or damage caused by use of our dialyzer after the expiry date mentioned on the label or packaging.

This device must be used on dialysis machines with an ultrafiltration controller or an accurate balancing system.¹

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



Materials²⁰
 Membrane²¹
 Housing²³
 Potting compound²⁵

Cellulose Triacetate²²
Polypropylene²⁴
Polyurethane²⁶

Sterilization²⁷
 Gamma Ray²⁸

Recommended connectors for blood ports Acc.to EN ISO 8637-1
 Recommended connectors for filtrate ports Acc.to EN ISO 8637-1
 Note: Operation of the dialyzer under clinical conditions may produce values different from those illustrated because of the variables involved in the clinical dialysis procedure, in the Cellulose Triacetate membrane, and in the manufacture of the device. Therefore, the values given are for approximate only.
 See in-vitro test conditions for explanatory materials relating to the test conditions from which the data were derived. Further information may be obtained on request.

EN FB-U PERFORMANCE DATA

Qb (mL/min)	05U	07U	09U	11U	13U	15U	17U	19U	21U
200	165	179	186	192	196	198	198	200	200
300	191	221	241	254	267	273	279	285	290
400	213	245	271	284	312	319	327	338	342
200	143	161	173	180	186	191	193	195	197
300	164	195	215	230	242	251	257	266	271
400	182	213	237	260	278	290	295	311	318
200	116	136	150	161	169	175	180	185	187
300	134	162	180	194	209	221	232	239	247
400	149	175	198	215	231	248	261	273	283
200	74	91	106	120	128	139	145	154	159
300	81	101	119	135	151	163	169	184	191
400	86	106	125	144	159	174	187	202	210
	15	19	21	24	29	33	37	41	45

Specifications ⁷	05U	07U	09U	11U	13U	15U	17U	19U	21U
Effective Surface Area ⁸ (m ²)	0,5	0,7	0,9	1,1	1,3	1,5	1,7	1,9	2,1
Priming Volume ⁹ (mL)	35	46	56	67	78	90	100	113	123
Effective Length ¹⁰ (mm)	162	173	188	202	213	224	233	242	251
Inner Diameter ¹¹ (µm)	200	200	200	200	200	200	200	200	200
Membrane Thickness ¹² (µm)	15	15	15	15	15	15	15	15	15
Maximum TMP ¹³ (mmHg)	500	500	500	500	500	500	500	500	500

Pressure Drops ¹⁴	Qb/Qd (mL/min)	200/500	200/500	200/500	200/500	200/500	200/500	200/500	200/500
Blood ¹⁵ /Dialysate ¹⁶ (mmHg)		67/8	58/12	52/13	49/11	46/11	44/11	41/9	39/10

In-Vitro Test Conditions¹⁷ (EN ISO 8637-1)
 Clearance : Qd 500mL/min, Qf 10mL/min
 KUF : Bovine Blood¹⁸ (Hct 32±2%, Protein¹⁹ 60g/L, 37°C), Qb 300mL/min

Materials²⁰
 Membrane²¹
 Housing²³
 Potting compound²⁵

Cellulose Triacetate²²
Polypropylene²⁴
Polyurethane²⁶

Sterilization²⁷
 Gamma Ray²⁸

Recommended connectors for blood ports Acc.to EN ISO 8637-1
 Recommended connectors for filtrate ports Acc.to EN ISO 8637-1
 Note: Operation of the dialyzer under clinical conditions may produce values different from those illustrated because of the variables involved in the clinical dialysis procedure, in the Cellulose Triacetate membrane, and in the manufacture of the device. Therefore, the values given are for approximate only.
 See in-vitro test conditions for explanatory materials relating to the test conditions from which the data were derived. Further information may be obtained on request.

EN ENGLISH	DE DEUTSCH	FR FRANÇAIS	ES ESPAÑOL	IT ITALIANO	NL NEDERLANDS	PT PORTUGUÊS	EL ΕΛΛΗΝΙΚΑ	SV SVENSKA	DA DANSK	NO NORSK
1 PERFORMANCE DATA	LEISTUNGSDATEN	DONNÉES DE PERFORMANCE	CARACTERÍSTICAS	CARATTERISTICHE	SPECIFICATIES	CARACTERÍSTICAS GERAIS	ΠΙΝΑΚΑΣ ΑΕΙΔΟΜΕΝΩΝ	PRESTANDA	PRESTANDA	VELESTEDATA
2 Clearance	Clearance/Werte	Clairance	Adiaramientos	Adiaramientos	Adiaramientos	Adiaramientos	Αδειασμα	Clairance	Clairance	Adiaramiento
3 Urea	Harnstoff	Urea	Urea	Urea	Urea	Urea	Ουρία	Urea	Urea	Urea
4 Creatinine	Kreatinin	Creatinine	Creatinina	Creatinina	Creatinina	Creatinina	Κρεατινίνη	Kreatinin	Kreatinin	Creatinina
5 Phosphate	Phosphat	Phosphate	Fosfato	Fosfato	Fosfato	Fosfato	Φωσφορικό οξύ	Fosfat	Fosfat	Fosfat
6 Vitamin B12	Vitamin B12	Vitamine B12	Vitamina B12	Vitamina B12	Vitamina B12	Vitamina B12	Βιταμίνη Β12	Vitamin B12	Vitamin B12	Vitamin B12
7 Specifications	Technische Daten	Spécifications	Características	Specifiche	Specificaties	Especificações	Χαρακτηριστικά	Specifikationer	Specifikationer	Specificasjoner
8 Effective Surface Area	Effektive Oberfläche	Surface Effective	Area Efectiva	Superficie Effettiva	Efficiënt Oppervlak	Area Efectiva	Επιφανειακό Εμβαδόν	Effektivi Yta	Effektivi Overfladeområde	Effektivi Flatenal
9 Priming Volume	Volumen	Volume d'Amorçage	Volumen de Cebado	Volumine	Volumine	Volume Estático	Χρόνος Εκκίνησης	Gründvolumen	Indvinkvolumen	Indvinkvolumen
10 Effective Length	Effektive Länge	Longueur Effective	Longitud Efectiva	Lunghezza Efficace	Effectieve Lengte	Comprimento Real	Εκτιμώμενο	Effektiv Längd	Effektiv Længde	Effektiv Lengde
11 Inner Diameter	Innendurchmesser	Diamètre Interne	Diametro Interior	Diametro Inter	Binnendiameter	Diametro Interno	Εσωτερικό Διάμετρο	Innendiameter	Indvendig Diameter	Indre Diameter
12 Membrane Thickness	Membrandicke	Epaisseur de Membrane	Espesor de la Membrana	Membrandikte	Membrandikte	Membrandikte	Πάχος μεμβράνης	Membran Tjkkelse	Membran Tjkkelse	Membrantjkkelse
13 Maximum TMP	Maximaler TMP	TMP Maximal	TMP Máximo	TMP Massimo	Maximale TMP	TMP Máximo	Μέγιστη TMP	Maximal TMP	Maximal TMP	Maximal TMP
14 Pressure Drops	Druckabfall	Perte de Charge	Caída de Presión	Calo di Pressione	Druckdaling	Diminuição de Pressão	Απώλεια πίεσης	Trykfall	Trykfall	Trykfall
15 Blood	Blut	Sang	Sangre	Sangue	Sang	Sangue	Αίμα	Blod	Blod	Blod
16 Dialysate	Dialysat	Dialysat	Dializato	Dializzato	Dialysaat	Solução de diálise	Αιματικό αμφοδιάλυτο	Dialysat	Dialysat	Dialysavevæske
17 In-Vitro Test Conditions	In-Vitro Testbedingungen	Conditions d'essai "In-Vitro"	Condiciones de Prueba "In-Vitro"	Condizioni di Prova "In-Vitro"	Testcondities In-Vitro	Condições de Teste "In-Vitro"	Συνθήκες δοκιμής in vitro	För In-Vitro-Test	För In-Vitro-Test	For In-Vitro-Test
18 Bovine Blood	Bovine Blut	Bovine Blood	Sangre Bovina	Sangre Bovino	Bosine Blood	Sangre Bovino	Βοσινικό αίμα	Storvet Blod	Storvet Blod	Storvet Blod
19 Protein	Protein	Protéines	Proteína	Proteine	Proteïne	Proteína	Πρωτεΐνη	Protein	Protein	Protein
20 Materials	Materialien	Matériaux	Materiales	Materiali	Materialen	Materialen	Υλικά	Material	Material	Material
21 Membrane	Membran	Membrane	Membrana	Membrana	Membran	Membrana	Μεμβράνη	Membran	Membran	Membran
22 Cellulose Triacetate	Celluloseacetat	Triacétate de cellulose	Triacetato de celulosa	Triacetato di cellulosa	Celluloseacetaat	Tricetato de celulose	Παράγωγο κυτταρίνης	Cellulosaacetat	Cellulosaacetat	Cellulosaacetat
23 Housing	Gehäuse	Carcasse	Carcasa	Corpo	Behuizing	Estrutura	Πλαστικό	Cellulosaacetat	Cellulosaacetat	Cellulosaacetat
24 Polypropylene	Polypolypropylen	Polypolypropylene	Polipropileno	Polipropilene	Polypolypropylen	Polipropileno	Πολυπροπυλένιο	Polypolypropylen	Polypolypropylen	Polypolypropylen
25 Potting compound	Verbindmittel	Résort de Fixation	Compuesto Adhesivo	Collante	Ingreetsmassa	Poliuretano	Πολυουρεθάνιο	Sammenstillning	Materiale sammensætning	Befestigelses sammensætning
26 Polyurethane	Polypolyurethan	Polyuréthane	Poliuretano	Poliuretano	Polyurethan	Poliuretano	Πολυουρεθάνιο	Polypolyuretan	Polypolyuretan	Polypolyuretan
27 Sterilization	Sterilisation	Sterilisation	Sterilización	Sterilizzazione	Sterilisatie	Esterilização	Αποστείρωση	Sterilisering	Sterilisering	Sterilisering
28 Gamma Ray	Gammastrahl	Rayons Gamma	Rayos Gamma	Raggi Gamma	Gammastraal	Raios Gamma	Ακτίνες γ	Gammestråling	Gammestråling	Gammestråling

EN SUREFLUX-E PERFORMANCE DATA

Qb (mL/min)	05E	07E	09E	11E	13E	15E	17E	19E	21E
200	154	172	182	191	196	198	199	200	200
300	187	215	234	247	263	270	279	284	289
400	207	242	268	275	292	306	314	323	328
200	130	152	165	176	182	188	192	193	195
300	153	181	200	214	229	237	244	251	257
400	166	201	224	239	253	268	278	287	296
200	110	128	142	149	161	166	173	179	181
300	119	148	165	178	192	201	213	222	230
400	124	154	176	199	216	229	244	256	264
200	59	74	87	97	106	115	123	130	135
300	63	81	96	101	116	125	133	143	151
400	67	86	100	108	123	132	143	154	166
	9	11	12	14	17	18	19	22	25

Specifications ⁷	05E	07E	09E	11E	13E	15E	17E	19E	21E
Effective Surface Area ⁸ (m ²)	0,5	0,7	0,9	1,1	1,3	1,5	1,7	1,9	2,1
Priming Volume ⁹ (mL)	35	45	56	68	78	90	101	113	122
Effective Length ¹⁰ (mm)	162	173	188	202	213	224	233	242	251
Inner Diameter ¹¹ (µm)	200	200	200	200	200	200	200	200	200
Membrane Thickness ¹² (µm)	15	15	15	15	15	15	15	15	15
Maximum TMP ¹³ (mmHg)	500	500	500	500	500	500	500	500	500

Pressure Drops ¹⁴	Qb/Qd (mL/min)	200/500	200/500	200/500	200/500	200/500	200/500	200/500	200/500	
Blood ¹⁵ /Dialysate ¹⁶ (mmHg)		66/11	55/13	52/13	49/12	46/12	44/11	42/11	41/10	40/10

In-Vitro Test Conditions¹⁷ (EN ISO 8637-1)
 Clearance : Qd 500mL/min, Qf 10mL/min
 KUF : Bovine Blood¹⁸ (Hct 32±2%, Protein¹⁹ 60g/L, 37°C), Qb 300mL/min

Materials²⁰
 Membrane²¹
 Housing²³
 Potting compound²⁵

Cellulose Triacetate²²
Polypropylene²⁴
Polyurethane²⁶

Sterilization²⁷
 Gamma Ray²⁸

Recommended connectors for blood ports Acc.to EN ISO 8637-1
 Recommended connectors for filtrate ports Acc.to EN ISO 8637-1
 Note: Operation of the dialyzer under clinical conditions may produce values different from those illustrated because of the variables involved in the clinical dialysis procedure, in the Cellulose Triacetate membrane, and in the manufacture of the device. Therefore, the values given are for approximate only.
 See in-vitro test conditions for explanatory materials relating to the test conditions from which the data were derived. Further information may be obtained on request.

EN SUREFLUX-G PERFORMANCE DATA

Clearance ² (mL/min)	Qb (mL/min)	11G	13G
---------------------------------	-------------	-----	-----

EN SUREFLUX™-L PERFORMANCE DATA¹

Qb (mL/min)	03L	05L	07L	09L	11L	13L	15L	17L	19L	21L
200	95	149	167	177	186	191	194	197	198	199
300	108	175	202	223	241	254	262	271	276	280
400	116	193	228	256	270	284	300	310	316	322
200	79	129	150	162	170	176	182	186	190	192
300	87	141	166	190	207	219	233	240	247	253
400	93	153	182	207	228	244	261	275	285	289
200	55	105	125	133	137	145	154	161	167	170
300	62	110	130	147	160	173	186	198	206	213
400	68	119	138	155	169	190	202	218	230	237
200	29	49	59	72	82	89	98	107	114	120
300	30	53	64	78	88	98	109	118	128	135
400	32	56	67	82	90	102	115	125	134	144
3	5	7	8	10	12	14	17	19	22	
03L	05L	07L	09L	11L	13L	15L	17L	19L	21L	
0,3	0,5	0,7	0,9	1,1	1,3	1,5	1,7	1,9	2,1	
22	35	45	56	68	78	90	100	113	123	
162	162	173	188	202	213	224	233	242	251	
200	200	200	200	200	200	200	200	200	200	
15	15	15	15	15	15	15	15	15	15	
500	500	500	500	500	500	500	500	500	500	
200/500	200/500	200/500	200/500	200/500	200/500	200/500	200/500	200/500	200/500	
111/14	67/10	57/11	53/13	50/11	49/11	46/9	45/11	42/11	40/9	

In-Vitro Test Conditions^{1,2} (EN ISO 8637-1)

Clearance : Qd 500mL/min, Qf10mL/min
 KUF : Bovine Blood^{1,8} (Hct 32±2%, Protein^{1,9} 60g/L, 37°C),
 Qb 300mL/min

Recommended connectors for blood ports Acclo EN ISO 8637-1
 Recommended connectors for filtrate ports Acclo EN ISO 8637-1

Note: Operation of the dialyzer under clinical conditions may produce values different from those illustrated because of the variables involved in the clinical dialysis procedure, in the Cellulose Triacetate membrane, and in the manufacture of the device. Therefore, the values given are for approximate only. See in-vitro test conditions for explanatory materials relating to the test conditions from which the data were derived. Further information may be obtained on request.

SK SLOVENČINA	SL SLOVENŠČINA	ET EESTI KEEL	AR العربية
1 UDAJE O VYKONE	PODATKI O UČINKOVITOSTI	ANALÜÜSI SELOOMUSTUS	بيانات الأداء
2 Clearance	Zračnost	Kliirens	التخفيض
3 Mochovina	Urea	Uurea	اليوريا
4 Kreatinin	Kreatiin	Kreatiiniin	الكرياتينين
5 Fosfaat	Fosfati	Fosfaat	فوسفات
6 Vitamiin B12	Vitamiin B12	Vitamiin B12	فيتامين ب 12
7 Spetsifikatsioon	Spetsifikatsioon	Spetsifikatsioonid	المواصفات
8 Efektivny povrch	Elektrivni del površine	Kasulik pindala	مساحة السطح الفعالة
9 Priming objem	Količina doziranja	Täitemah	حجم التهيئة
10 Efektivna alžka	Učinkovita dožina	Kasulik pikkus	الطول الفعّال
11 Vnitromy priemer	Notranji premer	Seesmine läbimõõt	القطر الداخلي
12 Hroúka membrány	Debelina membrane	Membrani paksus	سكس الغشاء
13 Maksimálne TMP	Maksimalni transmembranski tlak	Maksimaalne TMP	القصوى درجة حرارة
14 Poklesy tlaku	Padeoc tlaka	Rõhulangus	انخفاضات الضغط
15 Kivi	Kri	Veni	الورن
16 Dialyzaat	Dializati	Dialisaat	مغزى
17 Testovacie podmienky in vitro	Testni pogoji in vitro	In vitro analüüsi tingimused	شروط الاختبارات الخلية
18 Hovádzia krv	Vlovska kri	Veise veri	دم بقري
19 Bielekovina, protein	Protein	Proteiin	بروتين
20 Materijali	Materijali	Materjalid	مواد
21 Membrána	Membrana	Membran	غشاء
22 Triacetát celulozy	Triacetatna celuloza	Tselliuloosriasetaat	ناتلي السيتات السيلوزور
23 Kevl	Ložka	Kõpus	سيت
24 Polypropyleen	Polipropilen	Polüpropeen	مغزى بروبيلين
25 Puzdríaca hmota	Mehanski del filtra	Täiteaine	موكب بونفخ
26 Poluretaan	Põluretan	Põluretaan	بولوريان
27 Steariltsia	Steariltsia	Steariltsiamina	ستيرلستامينو
28 Gama luče	Gama žarki	Gammakiired	التعقيم

EN SUREFLUX™-N PERFORMANCE DATA¹

Qb (mL/min)	05N	07N	09N	11N	13N	15N	17N	19N	21N
200	134	152	166	175	181	186	191	193	194
300	155	179	201	216	230	239	248	256	261
400	165	195	220	241	256	268	279	291	303
200	105	122	137	149	158	165	169	175	178
300	119	139	158	175	188	196	207	215	222
400	129	156	179	198	213	223	236	249	261
200	68	84	97	109	120	129	136	143	147
300	75	90	107	120	133	140	151	160	168
400	82	97	114	129	141	152	162	174	183
200	32	41	49	55	62	68	73	80	86
300	35	43	51	58	66	71	77	84	92
400	37	45	54	61	68	74	81	88	95
3	4	5	6	7	8	9	10	11	
05N	07N	09N	11N	13N	15N	17N	19N	21N	
0,5	0,7	0,9	1,1	1,3	1,5	1,7	1,9	2,1	
35	46	57	66	77	88	100	113	124	
162	173	188	202	213	224	233	242	251	
200	200	200	200	200	200	200	200	200	
15	15	15	15	15	15	15	15	15	
500	500	500	500	500	500	500	500	500	
200/500	200/500	200/500	200/500	200/500	200/500	200/500	200/500	200/500	
64/9	55/11	50/13	48/10	45/11	43/9	41/10	38/9	36/9	

In-Vitro Test Conditions^{1,2} (EN ISO 8637-1)

Clearance : Qd 500mL/min, Qf 10mL/min
 KUF : Bovine Blood^{1,8} (Hct 32±2%, Protein^{1,9} 60g/L, 37°C),
 Qb 300mL/min

Recommended connectors for blood ports Acclo EN ISO 8637-1
 Recommended connectors for filtrate ports Acclo EN ISO 8637-1

Note: Operation of the dialyzer under clinical conditions may produce values different from those illustrated because of the variables involved in the clinical dialysis procedure, in the Cellulose Triacetate membrane, and in the manufacture of the device. Therefore, the values given are for approximate only. See in-vitro test conditions for explanatory materials relating to the test conditions from which the data were derived. Further information may be obtained on request.

NIPRO CORPORATION
 3-26, Senriokashinmachi, Settsu,
 Osaka, 566-8510, Japan

EC REP
NIPRO MEDICAL EUROPE
 Blokhuisstraat 42, 2800 Mechelen, BELGIUM

CH REP
Theramed AG Therapie- und Medizintechnik
 Sagihof 7, 6043 Adligenswil, Switzerland
NIPRO CORPORATION ODATE FACTORY
 8-7, Hanukiyachi, Niida, Odate-shi,
 Akita, 018-5794 JAPAN

[Distributed in Latin America by]
NIPRO MEDICAL CORPORATION
 3150 N.W. 107th Avenue, Miami, FL 33172,
 U.S.A.

[Authorized Representative in Malaysia]
NIPRO MALAYSIA SDN. BHD.
 B-3-2, The Ascent Paradigm, No.1 Jalan SS7/26A,
 Kelana Jaya, 47301 Petaling Jaya, Selangor,
 MALAYSIA
 +603-7622 1190, office@nipro.com.my
 Reg. No. for SUREFLUX-UX,
 E, G, L, and N: GC62461390017