SAFETY HEMODIALYSIS NEEDLE

〈Instructions for Use〉

Prior to use, please read the Instructions for Use carefully.

I. Indications

Safety Hemodialysis Needle is used in the dialysis treatment of patients with acute or chronic renal failure. Blood access will be secured by a healthcare professional.

II. Contraindications

- Do not re-use.
- Do not push forward the inner needle inside of the plastic cannula once the inner needle has been partially pulled out. It may damage the plastic cannula or puncture it. Do not move the inner needle forward or backward. [The plastic cannula or hemostasis valve may break and blood may leak (See Fig. 1).]

III. Operating Principle



The stopper of the safety mechanism gets stuck with the protrusive area of the inner needle when the inner needle is retracted from the plastic cannula after having punctured the patient. At the same time, the magnet inside the safety protector shifts to occlude the small hole of the safety protector and protects the bevel.

IV. Clinical Benefits

Dialysis can be performed by connecting the device with a dialysis machine, blood circuit, and hemodialyzer to maintain the life of renal patients.

V. Warnings

1. Connect the luer connectors of the plastic cannula and the other device straight. When they are connected askew, the connectors may be disconnected during use due to a weak engagement force and hazardous situations such as a blood leak may produce symptoms such as a decrease in blood pressure in patients.
2. Medication such as lidocaine air spray or alcohol on the devices may cause hazardous situations such as a blood leak due to a crack, which may produce symptoms such as a decrease in blood pressure in patients.



VI. Cautions

1. Do not use the device if the package is broken, or an abnormality like a broken device, etc. is found.
2. Re-use of the device may lead to infection by contamination.
3. Use immediately after opening, be careful to prevent infection after use and dispose of the device in a safe method.
4. Pull the inner needle out straight and do not bend the plastic cannula. Do not touch the safety mechanism. If the safety mechanism is touched or the inner needle is pulled out askew, the safety mechanism may not function. In case that the safety mechanism is detached from the luer cap before pulling the inner needle, dispose of the inner needle in a safe method.
5. Do not touch the safety mechanism attached to the inner needle tip that was drawn from the plastic cannula; dispose of it promptly with the inner needle. Pulling or scratching the safety mechanism may cause it to be separated or blood may scatter. Applying excessive pressure may cause the needle tip to expose itself from the safety mechanism.
6. Do not use a bent inner needle. If the inner needle has been pulled out with the inner needle bent, the safety mechanism against needle stick injuries may not work.
7. Use caution on the crack of the connector when injecting fat emulsion, pharmaceutical products containing fat emulsion, oily ingredients such as castor oil, solubilization agents such as surface active agent or alcohol, or using disinfectants containing alcohol. [Medication may cause cracks in the connector, blood and medication leaks or air inclusion. Especially note that if the required injection volume is not obtained due to these accidents in the injection of full anesthesia, vasopressor, anticancer drugs and immunosuppressive drugs, this may cause serious effects on the patients. Re-fastening at the time of line exchange, excessive fastening or additional fastening may cause cracks.]
8. Magnetic material is used in the safety mechanism. Thus, moving the device closer to precision equipments containing magnetic material of a heart pacer, or a magnetic recording medium, etc., the device may malfunction. Take caution if using those equipments.
9. Do not touch the needle tip of the device directly.
10. Use caution not to adhere medication or blood to the connecting section. [Loosening of the connecting section may occur.]
11. Do not use the device if abnormalities such as blood leak from the connecting section or air inclusion, etc. are found.
12. Do not apply excessive force, pressure, etc. to the device. [The device may be damaged or deformed.]
13. Check for looseness or a liquid leak at the connecting section before use. During use, check periodically for device damage, loosening at the connecting section, blood leaks, air inclusion and clogging.
14. Connections and disconnections between the device and other devices must be performed quickly. [If other devices are connected or disconnected slowly, blood may leak from the hemostasis valve. Also, if they are removed forcefully, the hemostasis valve may leak blood until it is resealed.]
15. If the hemostasis valve remains open due to friction or blood clots, clamp the tube while avoiding contact with the leaking blood.
16. After ensuring that the flow guide has passed completely through the hemostasis valve and the flow is sufficient, start the treatment by extracorporeal circulation.
17. Do not fasten the connector with an excessive force. [Do not separate or damage the connector, or liquid leaking from the connecting section and air inclusion may occur.]
18. Dispose of the device in an approved biohazard container as per facility protocol.
19. Any incidents shall be reported to the manufacturer and the competent authority of your State.
20. Do not store at extreme temperature and humidity.
21. If the device is exposed to abnormal conditions (for example, high temperature and humidity) or if it is unintentionally opened before use, do not use the device.

VII. Instructions for Use

1. Conduct all operation aseptically.
2. Disinfect the puncture site. <Note> Shave hair before disinfection if necessary.
3. Pull the device out straight from the package. <Caution> Open the package just before use. When removing the device from the package, pull the device out straight. If the device is pulled out askew, the inner needle tip may be damaged.
4. Check if the plug is firmly inserted into the inner needle hub. (See Fig. 3)
5. Re-tighten the luer cap gently. Check if the luer cap is secured tightly enough not to leak blood during puncturing. (See Fig. 4)
6. Check if the inner needle and plastic cannula are set at the correct position. <Caution> Do not use the plastic cannula if it protrudes from the inner needle.
7. Hold the inner needle hub and puncture the site. If blood inflows into the flash back section, it is verified that the plastic cannula has been inserted inside a blood vessel (vascular access secured). <Caution> Hold the inner needle hub when puncturing. If you puncture while holding other sections, or without retaining the inner needle, the inner needle may move during puncture, and thus it may not be so easy to puncture.
8. Hold the plastic cannula in one hand, and pull out the inner needle with your other hand. <Caution> When pulling out the inner needle, keep the plastic cannula straight. If it is found that vascular access has not been secured while pulling out the inner needle, pull out both the inner needle and plastic cannula and dispose of them. <Caution> Do not strongly nip the plastic cannula with sharp instruments such as forceps or your nails. The plastic cannula may be broken or damaged.
9. Check if the safety mechanism functions and inner needle tip is securely protected. (When the operation of safety mechanism completes, the inner needle does not return to original position.) <Caution> When removing the inner needle without protecting the inner needle tip, be careful of wrong puncture.
10. Insert the plastic cannula inside the blood vessel up to the required depth.
11. Remove the luer cap and connect the luer connectors of the plastic cannula and the other device straight. <Caution> Make sure that the plastic cannula and the other device are connected securely with luer connectors. Verify that the screw of the lock mechanism is securely fastened.
12. During dialysis, fix the luer connectors with tape and/or cover the puncture site with gauze if necessary.

VIII. Connection with other devices

The device connects to a blood circuit or other circuit for dialysis, by using a connector with a luer taper.

IX. Guarantee

1. Our Safety Hemodialysis Needle is manufactured under strict quality control conditions and the quality is assured. We shall not be responsible for injury caused to a patient or to any person or damage to any object that is attributed to transport, storage, and usage in your institution.
2. If a patient or any person is injured or any object is damaged when using our Safety Hemodialysis Needle, we shall not be responsible for the injury or damage unless we are clearly identified as being at fault.
3. If a patient or any person is injured or any object is damaged when re-using our Safety Hemodialysis Needle, we shall not be responsible for the injury or damage of any nature.
4. We shall not be responsible for any injury or damage caused if our Safety Hemodialysis Needle is used after the expiry date mentioned on the packages.

[Symbols used for labeling]

STERILE R	Sterilized using irradiation
	Contains hazardous substances

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