SOLACEA™

HEMODIALYZER

Instructions for Use

Indications: Hemodialysis or Hemodiafiltration with SOLACEA-H is indicated for patients having acute or chronic renal failure when dialysis is prescribed by the physician.

Contraindications: Do not re-use. Patients indicating allergic reactions to cellulose triacetate membranes should not be dialyzed with this product.

Clinical Benefits: Dialysis treatment can be performed by connecting the device to a dialysis machine, a blood line, and an AVF needle, and waste and unnecessary water is removed from the blood to support the life of renal failure patients.

Warning

- Use this product according to instructions of a physician who is well familiar with the patient's condition.
- Follow these instructions and those of the dialysis machine supplier.
- Do not use for any other purposes than dialysis.
- If any abnormalities such as foam generation or mixture, blood leakage, blood coagulation and hemolysis occurred during the use of this product, take appropriate measures according to a physician's instructions.
- In case drugs including an anticoagulant are administered before or during use of this product, follow a physician's instructions about the administration and dose and the administration time of the drugs.
- Do not re-use this product since this is a single-use product.
- The foreseeable risks in association with re-use of the product are:
 - Infection by contamination,
 - Deterioration of solute removal performance and ultrafiltration performance,
 - Exposure of patients and/or technicians to residual medicinal agents such as disinfectant used for product re-use, and/or adverse effects of residual medicinal agents on them, and
 - Damage of hollow fiber and/or leakage.
- Do not expose this product to chemical solvents, such as bleach and alcohols. Residual disinfectant in this product may cause adverse patient reactions.
- If the patient exhibits any abnormal symptoms such as discomfort, pruritus, urticaria, peripheral and facial edema, respiratory arrest, facial flush, erythema, asthmatic reaction, hypertension, hypotension and/or cardiac arrhythmia during the use of this product, take appropriate measures according to a physician's instructions.
- Commonly seen side effects (hypotension, hypertension, headache and nausea which are sometimes with hypovolemia or hypervolemia) can be minimized by careful management of the patient fluid and electrolyte balance, as well as the dialysis condition (blood flow rate and ultrafiltration rate).
- During dialysis, constantly monitor the patients who ;
 - (1) have a history of hypotension with therapy.
 - (2) have inflammatory reaction, allergic reaction, hypersensitivity, or increase in the immunity by infections.
 - (3) take hypotensive drugs such as inhibitor of angiotensin converting enzyme and calcium antagonist.
 - (4) use this product for the first time.

Caution: Caution should be employed against excessive water removal. Use of accurate UF control system is required. Do not use on non-de-aerated dialysis fluid delivery systems. Confirm that the dialysate does not contain pyrogens in order to prevent transfer of pyrogens from dialysate to blood.

Any incidents shall be reported to the manufacturer and the competent authority of your State.

Caution before use

- Do not use if package is damaged.
- Do not use the device if the protective end caps are not in place.
- Unpack immediately before use.
- Avoid air mix-in and contamination during rinsing / priming operations.
- Start dialysis immediately after rinsing / priming operations.
- Rinsing / priming should be carried out under the following conditions according to this "Instructions for Use":
 - Blood side : Rinse with physiological saline (not less than 1 000 mL) at a flow rate of 100 mL/min.
 - Dialysate side : Verify conductivity and temperature, and rinse with dialysate at a flow rate of 500 mL/min for about 5 minutes.
- Check the integrity of the blood line and dialyzer.
- Administration of Heparin
Systemic or regionalized heparinization may need to be administered based on instructions from the attending physician.

Caution in use

- Continuously monitor the pressure in the blood line and check for blood leakage during dialysis.
- Carefully avoid contamination during blood sampling and blood reinfusion.
- Set TMP alarm (max. 500 mmHg).
- Avoid air during blood reinfusion to minimize the risk of air embolisms.
- Do not apply excessive pressure to the blood line, the dialyzer and their connections.

Caution after use

- Single use only. Dispose of the dialyzer immediately after use.
- Dispose of the used blood lines and dialyzer by any means suitable for avoiding contamination.
- Dispose of the device in an approved biohazard container as per facility protocol.

Caution for storage

- If the device is exposed to abnormal conditions (for example, high temperature and humidity) or if it is unintentionally opened before use, do not use the device.
- Handle carefully in low temperature environment. (Especially below freezing, it may result in breakage of the hollow fiber.)

Instructions for Use

I. Rinsing / priming

Follow instructions on the machine's operator manual.

«In case of priming with a physiological saline bag»

- (1) Take the dialyzer out of the package and set it to the holder so that the venous side is directed upward. (Fig. 1)
- (2) Connect the arterial and venous dialyzer connector to dialyzer. Connect the arterial patient line to physiological saline bag (Fig. 2). Start arterial blood pump at a flow rate of 100 mL/min (not less than 1 000 mL) (Fig. 3).
- (3) Stop the blood pump; rotate the dialyzer for 180 degree. Put the dialysate connectors to the dialyzer (dialysate inlet at the venous blood side, dialysate outlet at the arterial blood side) (Fig. 4). Make sure of degassing the dialysate part of the dialyzer. Start running through dialysate at a flow rate of about 500 mL/min.
- (4) Restart the blood pump. Make sure that the blood compartment is free of air-bubbles and filled with physiological saline. The preparations for dialysis is complete (Fig. 5).
«In case of online priming»
Follow instructions on the machine's operator manual.

[Leakage test] In case of physiological saline priming

It is recommended to perform the following operations before connecting the dialysate lines to the dialyzer.

- (1) Fully prime the arterial and venous lines and the dialyzer with physiological saline by operating the blood pump ; then stop the pump operation.
- (2) Clamp the arterial line near the dialyzer and the distal end of the venous line with clamps.

- (3) Place the clamped distal end about 1 m below the dialyzer and remove the clamps. (This results in application of a negative pressure of approximately 70 mmHg to the blood compartment of the dialyzer.)
- (4) Examine whether or not continuous bubble formation is observed in the venous header to check for leakage from the dialyzer, if observed, replace the dialyzer with a new one.

II. Start of dialysis

Follow instructions on the machine's operator manual.

- (1) Prepare the blood access site and connect to the arterial line. Remove the clamps from the arterial and venous lines. While running dialysate at a flow rate of approximately 500 mL/min, operate the blood pump at a flow rate of approximately 50 mL/min.
- (2) Confirm that no air bubbles remain in venous header or venous blood line.
- (3) Fully prime the arterial and venous lines including the dialyzer with blood by operating the blood pump, and then stop the pump operation. Clamp the distal end of the venous line with clamps.
- (4) Prepare the blood return site and connect to the venous line. After confirming that there are no bubbles in the line, remove the clamps from the line. After checking that there are no clamps on the lines and no line folding, operate the blood pump at a low flow rate.
Check the integrity of the connections.
- (5) After confirming that there are no bubbles in the arterial and venous headers, turn the dialyzer 180° to allow removal of bubbles from dialysate. If bubbles are detected in the venous header before the turning, run blood at a prescribed flow rate for 5 to 10 minutes with the venous side kept upward.

III. Operations during dialysis

- (1) If stopping the blood pump is required during dialysis due to insufficient blood flow, lower the dialysate pressure to approximately 0 mmHg. (This is to avoid blood coagulation due to dehydration.)
- (2) Set the UF rate carefully to avoid excessive water removal according to the patients' needs. Reduce the blood flow rate if disequilibrium syndrome is suspected.
- (3) If blood leakage is suspected, judge by testing dialysate sampled from the dialysate outlet port of the dialyzer using occult blood reaction test paper. If a leakage is detected, reduce the UF rate to minimum rate according to the institutional protocol, stop the dialysate supply and blood reinfusion, then replace the dialyzer with a new one.

IV. Dialysis termination and blood reinfusion

- (1) Stop the blood pump, clamp the arterial line and remove the line from the arterial blood access site ; then connect the line to the physiological saline vial for blood reinfusion.
- (2) Unclamp the arterial line and run physiological saline to rinse out the blood from the arterial and venous lines and the dialyzer.
- (3) After the blood reinfusion, discard the arterial and venous lines and the dialyzer. Do not reuse them.

Performance and Specification

The performance of the hemodialyzer varies with types.

Refer to respective catalogues and performance data sheets.

Connection with other devices

The device contacts and is connected with the blood line at the header and the dialysis machine through the coupler port of the housing.

Guarantee

- Non-pyrogenic, sterile fluid path that has been sterilized using irradiation.
- This dialyzer is manufactured under strict quality control and the quality is assured. If the dialyzer is defective (broken package, damaged dialyzer), it shall be replaced with a new one at our cost upon return of the broken package or damaged dialyzer. We will not be responsible, however, for the injury of a patient or any person or the damage to any object that is attributed to transport, storage and operation in your institution.
- If a patient or any person is injured or any object is damaged by use of our dialyzer, we will not be responsible for the injury or damage unless we are clearly identified as being at fault.
- If a patient or any person is injured or any object is damaged by reuse of our dialyzer, we will not be responsible for the injury or damage of any nature.
- We will not be responsible for any injury or damage caused by use of our dialyzer after the expiry date mentioned on the label or packaging.

This device must be used on dialysis machines with an ultrafiltration controller or an accurate balancing system.¹

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




ATA™
Cellulose Triacetate
Polypropylene²⁶

Recommended connectors for blood ports Acc, to EN ISO 8637-1
Recommended connectors for dialysate ports Acc, to EN ISO 8637-1

Note: Operation of the dialyzer under clinical conditions may produce values different from those illustrated because of the variables involved in the clinical dialysis procedure, in the ATA Cellulose Triacetate membrane, and in the manufacture of the device. Therefore, the values given are for approximate only. See in-vitro test conditions for explanatory materials relating to the test conditions from which the data were derived. Further information may be obtained on request.

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NIPRO CORPORATION
3-26, Senriokashinmachi, Settsu,
Osaka, 566-0810, Japan

EC REP
NIPRO MEDICAL EUROPE
Blaikins Road, 42, 2800 Mechelen, BELGIUM

CH REP
Theramed AG Therapie- und Medizintechnik
Siglist 7, 6043 Adligenswil, Switzerland

(Manufacturing facility)
NIPRO CORPORATION ODATE FACTORY
8-7, Hanukiyachi, Nida, Odate-shi, Akita,
018-0794 JAPAN

(Distributed in Latin America by)
NIPRO MEDICAL CORPORATION
3150 N.W. 107th Avenue, Miami,
FL 33172, U.S.A.

(Sponsor in Australia)
NIPRO AUSTRALIA PTY LTD
Suite 202, Level 2, 657 Pacific Highway,
St Leonards, NSW 2069, AUSTRALIA

(Distributed in Brazil by)
NIPRO MEDICAL LTD.
Avenida Nipro IV-651, Região Norte
CEP 18.087-127 Sorocaba SP BRASIL

DCE-AK-SOLACEAMDR-address=00A

DCE-AK-SOLACEAMDR-00E (Creation date: 2024-11-05)