

Arterial cannule



Special features

- Arterial cannule has no step between the tip and body to allow smooth blood sending.
- Kink-resistant wire-reinforced body with depth marks.
- Thin tip wall contributes to maintain low pressure drop of blood.

Arterial cannule(45 degrees tip)

Specifications

| Code | Diameter [Fr] | Length [mm] | Connector [inch] |
|-------------|---------------|-------------|------------------|
| CAPA45+18CT | 18 | 341 | |
| CAPA45+20CT | 20 | | 3/8 |
| CAPA45+22CT | 22 | 343 | |
| CAPA45+24CT | 24 | | |