



Arterial Venous Fistula Needle (AVF)

High Quality Performance with Extra Thin Wall Needle

Specifications

- Wing : Turnable Wing or Fixed Wing
- Cannula : Extra thin wall with / without back-Eye

- Tube Length : 300mm
150mm
- Sterilization : ETO / Gamma

Connector & Clamp

Color	Code
White	W
Red	R
Blue	B



Product Range

Gauge	Needle Length	Tube Length	Wing Color		Remarks
			NIPRO	ISO	
14G	4/5" (20mm)	150mm or 300mm	Violet	Pale-Green	with Back Eye, Clamp, Turnable or Fixed Wing available upon your request
	1" (25mm)				
	1-1/4" (32mm)				
15G	4/5" (20mm)	150mm or 300mm	Ivory / Blue*	Blue-Grey	
	1" (25mm)				
	1-1/4" (32mm)				
16G	4/5" (20mm)	150mm or 300mm	Light-Green	White	
	1" (25mm)				
	1-1/4" (32mm)				
17G	4/5" (20mm)	150mm or 300mm	Orange	Red-Violet	
	1" (25mm)				
	1-1/4" (32mm)				

Package : 500 Pcs / Carton for 270mm Tube Length. 1000 Pcs / Carton for 120mm Tube Length.

* for USA market

* Single needle available with needle length 1" (25mm)