


BioHole™ DIALYSIS CATH



Special Facet Configuration Designed for Buttonhole Technique



Specifications

- Sterilization : Gamma
- Flow Rate : 16G : -185ml / min
17G : -153ml / min

CE 0123

Product Range

Package : 1,000 Pcs / Carton

Catheter		Needle		Color Code
Gauge (mm)	Length	Gauge (mm)		
16G (1.6mm)	25mm	18G (1.3mm)		Yellow-Green
	32mm			
	38mm			
17G (1.4mm)	25mm	19G (1.1mm)		Orange
	32mm			
	38mm			