



Dull AVF (BioHole Needle)

Special Facet Configuration Designed for Buttonhole Technique

Specifications

- Wing : Turnable Wing or Fixed Wing
- Cannula : Enables less damage on vascular access to buttonhole patients
- Tube Length : 400mm
300mm
150mm
- Sterilization : ETO / Gamma



Product Range

Gauge	Needle Length	Tube Length	Wing Color		Remarks
			NIPRO	ISO	
14G	4/5" (20mm) 1" (25mm)	150mm, 300mm or 400mm	Violet	Pale-Green	with Back Eye, Clamp, Turnable or Fixed Wing available upon your request
15G	4/5" (20mm) 1" (25mm)	150mm, 300mm or 400mm	Ivory / Blue*	Blue-Grey	
16G	4/5" (20mm) 1" (25mm)	150mm, 300mm or 400mm	Light-Green	White	
17G	4/5" (20mm) 1" (25mm)	150mm, 300mm or 400mm	Orange	Red-Violet	

Package : 500 Pcs / Carton for 300mm and 400mm tube length. 1000Pcs / Carton for 150mm tube length.

* for USA market