

FA-40



Arterial Filter

High Performance of Eliminating Particles and Air Bubble from the Blood during Extracorporeal Circulation

Special Features

- Sterilization : ETO

Specifications

| | |
|------------------------------|---------------------|
| Priming | 220 mL |
| Retention pores | 40 μ m |
| Effective filtering surface | 418 cm ² |
| Inlet and outlet | 9.53 mm (3 / 8") |
| Maximum flow | 7L / min |
| Connector for 3-way stopcock | Luer lock type |

Materials

| | |
|------------|-------------------------------|
| Filter | Monofilament polyester fabric |
| Body | Polymethylmetacrylate |
| Lines | PVC |
| Connectors | Polymethylmetacrylate |
| Casing | Polyurethane |