



# Safetouch™ Dialysis Cath Plus

Preventing Needle Stick Injury with Reliable Safety Function Specifications  
Preventing Blood Spreading with Auto Priming Hemostatic Valve

- Sterilization : Gamma
- Passive Safety Technology
- Flow Rate : 16G×38mm : 202.9mL / min



## Product Range

Package : 500Pcs / Carton

Catheter		Needle	hub color
Gauge(mm)	Length	Gauge(mm)	
15G(1.8mm)	25mm	17G(1.4mm)	blue gray
	32mm		
	38mm		
16G(1.6mm)	25mm	18G(1.2mm)	white
	32mm		
	38mm		
17G(1.4mm)	25mm	19G(1.1mm)	redviolet
	32mm		
	38mm		