



SAFETOUCH™ TULIP AVF

Safety Arterial Venous Fistula needle set with Anti-Stick needle Protector Features of the clear, Patented TULIP safety needle protector

Specifications

- Easy of use during patient connection and needle placement
- Activation using a one-handed technique and meet OSHA safety requirement
- Audible click telling protector is locked and will stay locked
- Cannula : Extra thin wall with/without back-eye
- Sterilization : ETO / Gamma



Product Range

Gauge	Needle Length	Tube Length	Wing Color		Remarks
			NIPRO	ISO	
14G	1" (25mm)	150mm or 300mm	Violet	Pale-Green	with Back Eye, Clamp, Turnable or Fixed Wing available upon your request
	1-1/4" (32mm)				
15G	1" (25mm)	150mm or 300mm	Ivory / Blue*	Blue-Grey	
	1-1/4" (32mm)				
16G	1" (25mm)	150mm or 300mm	Light-Green	White	
	1-1/4" (32mm)				
17G	1" (25mm)	150mm or 300mm	Orange	Red-Violet	
	1-1/4" (32mm)				

Package : 500 pcs of single pack per box. 120 pairs plus 10pcs of single pack per box.

* for USA market