



## UNIQUE

### Cardioplegia System

#### Special Features

- Flexibility and safety for a wide range of use
- Effective performance in a compact low priming design
- Pressure manometer and temperature probe port
- High heat exchanger performance
- Sterilization : ETO

#### Specifications

Heat exchanger effective area	0.02 m <sup>2</sup>
Heat Exchanger priming volume	30 mL
Equipment priming volume	105 mL
System total volume	135 mL
Crystalloid solution filter pore diameter	0.8 microns
Filtering area crystalloid solution filter	14.5 cm <sup>2</sup>

#### Materials

Heat exchanger	Anodized Aluminium
Central Valve	Polyethylene