

BRIZIO™ Arterial filter

Arterial Filter

High performance of eliminating particles and air bubbles from the blood during extracorporeal circulation

Special Features

- Sterilization : ETO



Specifications

Priming	96 mL
Filter opening size	40 µm
Filter opening area	35%
Effective filtering surface	450 cm ²
Blood inlet and outlet	3/8 inch
Maximum blood flow rate	7 L/min
Connector for 3-way stopcock	Luer lock type

Materials

Filter	Monofilament polyester fabric
Body	Polycarbonate
Lines	Polyvinyl chloride
Connectors	Polymethyl methacrylate
Casing	Polyurethane

Model code

Code	Description
CDD38	Arterial filter with bypass line