

# SURDIAL™ 55Plus

Compact and simple designed to allow easy settings and operations

## Special Features

- Compact body size, 280mm width
- 180 degrees swivable touch panel display
- Single needle system with one clamp
- UF profile and Conductivity profile
- Detection of micro air bubbles

CE<sub>0123</sub>

## Specifications

General data	
Dimension	280(W)×420(D)×1365(H)
Weight	Approx. 72kg
Power Supply	AC230V±10%, 50/60Hz AC110V±10%, 50/60Hz
Power consumption	Max. 2000VA
Water Supply	Temperature : 15 to 30°C 5 to 30°C (with heat exchanger, option) Pressure : 0.05 to 0.74 MPa Flow rate : 900mL or more
Extracorporeal circuit	
Arterial pressure monitoring	Measurement range : -500 to 500 mmHg Accuracy : ±10mmHg
Venous pressure monitoring	Measurement range : -500 to 500 mmHg Accuracy : ±10mmHg
Air Bubble Detector	Ultrasonic Sensor Two conditions create an air bubble alarm: any single bubble greater than 10 µL in size and when bubbles of 0.3 µL or larger, accumulate over a 1-minute period into a volume of 1, 10, 50 or 300 µL.
Arterial blood pump	Flow range : 0,10 to 600mL/min Accuracy : ±10%
Heparin pump	Delivery range : 0.0 to 20.0 mL/hr Bolus range : 0.1 to 9.9 mL Accuracy : ±10% Syringe size : 10,20,30mL
Hydraulic circuit	
Ultrafiltration	Volumetric control UF rate:0,0.10 to 5.00L/hr Accuracy: ±30g/hr
Dialysate flow rate	300 to 800 mL/min
Dialysate temperature	34.0 to 39.0°C

Concentration range	Bicarbonate dialysis Water +B : 5 to 65 mmol/L Water +A+B : 130 to 160 mmol/L (Na <sup>+</sup> electrolyte concentration) Stability : 0.3mS/cm or less Acetate dialysis Water +A : 130 to 160 mmol/L (Na <sup>+</sup> electrolyte concentration) Stability : 0.3mS/cm or less
Blood leak detector	Optical
Rinse Program	Water rinse Chemical disinfection Chemical decalcification Hot rinse Hot citric disinfection Tepid rinse
Function	
Acetate and Bicarbonate dialysis Isolated ultra filtration Single needle with single blood pump UF profile Concentration profile Dose finder (Kt/V calculator) 10.4" color touch display	
Option	
Blood pressure monitoring (BPM) LAN connector Front hook for solution bag CF flushing unit Dialysate solution tank tray Acid rinsing port unit Potential equalization conductor * Heat exchanger unit(Factory Option) * Bicarbonate powder cartridge system * Central system for A solution (Factory Option)	

\*: Factory installed option